## OPHTHALIGHT OPHTHALMIC DEVICES



ML CI is an Excellent IOL delivery system that utilizes the innovative treatment technology specifically designed and developed to enable a safer, simpler and more effective intraocular lens (IOL) delivery, which protects the IOL from damage during cataract surgery. The ML CI IOL Delivery System, packaged in sets of One, is available in a range of incision sizes. Each set contains a sterile single-use injector and a treated cartridge.

Model reference	External diameter	Clear corneal incision using wound assisted	Intra cameral	Recommended IOL
		technique (dock-in)	incision	powers
ML CI 28	1.85	2.2	2.6	less than 24.0 D
ML CI 32	2.10	2.4	3.0	Higher than 24.0 D





No.1468, Derakhshan Ave., Shahrak Gharb, Tehran, Iran Tel.: 0098 21 88092641 www.mehrdavar.com